



# Certificate of Analysis

EXCELLENT **CUSTOMER CARE** IS AND CONTINUES TO BE THE CORE PHILOSOPHY OF INNOVATIVE RESEARCH. OUR **COMMITMENT** TO OUR CUSTOMERS IS REFLECTED IN HOW WE HAVE BUILT OUR NETWORK, THE PRODUCT RANGE WE OFFER, THE SERVICES WE PROVIDE AND HOW QUICKLY WE RESPOND TO CUSTOMERS' FEEDBACK.

<b>Product:</b>	Beagle Complement Serum
<b>Lot Number:</b>	Sample
<b>Number of Units:</b>	Sample
<b>Volume of Units:</b>	Sample
<b>Description:</b>	Pooled low hemolysis serum recovered from non-medicated, non-immunized, spontaneously clotted whole blood. Centrifugation was used to separate the serum from the clot. Material flash frozen to preserve complement activity. Complement Activity: 1:80 microtiter plates
<b>Storage:</b>	-20°
<b>Appearance:</b>	Amber/ Red colored liquid
<b>Expiration:</b>	Animal Plasma and serum is stable up to three years from the date of receipt when stored at - 20° or colder and not subject to multiple freeze thaw cycles.

This product is for research use only. This material should be handled as if capable of transmitting infectious agents. Please use universal precautions.

All materials are obtained from FDA licensed facilities (USA). Store at -20° C or colder

The reagents that you purchase from Innovative Research, Inc are clearly labeled "For Research Use Only". To comply with U.S. Food and Drug Administration (FDA) Regulations, these products are Not for use in Clinical, Diagnostic or Therapeutic Procedures. As your supplier, we advise our customers and monitor the use of these products to ensure that they are used for research purposes only. If you have any questions, do not hesitate to contact us. Thank you for your interest in Innovative Research Inc