



## Certificate of Analysis

PRODUCT:	Ultra Pure Rabbit Gamma Globulin
LOT NUMBER:	Sample
NUMBER OF UNITS:	Sample
VOLUME OF UNITS:	grams
PROCEDURE:	Ultrapure rabbit gamma globulin is purified from normal rabbit serum (either sex, healthy animals) using different fractionation procedures and extensive dialysis to remove insoluble lipoproteins. The product is filter through a 0.2 micron filter and frozen. Following immunoelectrophoresis versus anti-whole serum, the product displays precipitin bands only in the gamma region.
STORAGE:	2-8 C
ASSAY:	<p>This lot contains &gt;98% albumin free, which is less than 2% impurities as determined be serum protein electrophoresis.</p> <p>Protein by Weight: O.D 280 nm E 1% = 0.66 = 90%. Residual salts, lipids, etc. comprise the rest of the weight.</p>
DATE:	Sample

This product is for research use only. This material should be handled as if capable of transmitting infectious agents. Please use universal precautions.

The reagents that you purchase from Innovative Research, Inc are clearly labeled "For Research Use Only". To comply with U.S. Food and Drug Administration (FDA) Regulations, these products are Not for use in Clinical, Diagnostic or Therapeutic Procedures. As your supplier, we advise



our customers and monitor the use of these products to ensure that they are used for research purposes only. If you have any questions, do not hesitate to contact us. Thank you for your interest in Innovative Research Inc